SHERAPRESS

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en INSTRUCTIONS FOR USE

Denture base acrylic, cold cure polymerisate based on methyl methacrylate in powder and liquid form. Free of cadmium, phthalates and BPA.

1. Indications

For the completion of full dentures using the casting technique, completion of cast partial dentures, total and partial relinings (indirect method), supplementation and repair works. Classification as per DIN EN ISO 20795-1 type 2 class 1 and as per Medical Devices Directive 93/42/EEC Annex IX class Ia for removable and fixed prostheses. Only for intended use by qualified personnel (dental laboratory/dental practice).

2. Contraindication

Patient contact with the unpolymerised acrylic or the individual components is contraindicated. With a confirmed allergy to one or more of the components of the medical device, refrain from using the device.

3. Properties

Comfortable range of processing options, very good flowability, optimal mucosal compatibility thanks to low residual monomer content, more homogeneous surface, very low susceptibility to plaque development, easy to finish and polish, absolute colour stability thanks to catalyst system with no tertiary amines, tested and certified biocompatibility.

4. Safety instructions

Contains methyl methacrylate; 1,4-butanediol dimethacrylate; dibenzoyl peroxide; liquid and vapour highly flammable, burns possible. Causes skin irritation. May cause an allergic skin reaction. May cause respiratory irritation. Avoid inhalation of dust / fume / gas / mist / vapour / aerosol. Wear protective gloves / protective clothing / eye protection / face protection. In case of contact with skin: wash with plenty of soap and water. Observe further instructions in the safety data sheet! Finished medical product: Flavour irritation or injuries due to breakage of the product possible. Contact with alcohol, coffee, tobacco smoke or certain medications can lead to discolouration.

*The times stated refer to a material and room temperature of 23°C. Higher temperatures shorten these times while lower temperatures extend them.

**The time stated of approx. 7 minutes requires a previous processing time of approx. 8 minutes. The time in the pressure vessel should last up to the 15 minutes after the mixing.

5. Processing instructions

Preparation

Fix teeth set up in wax using a plaster or silicone matrix. Then carefully clean the matrix, model framework and the teeth by boiling out any residual wax (e.g. by using SHERAWAX-EX) and then steam off. The plaster model is wetted until no more bubbles rise up so that during pressure polymerisation and when allowing the acrylic to flow or when pressing, no air is forced from the plaster into the acrylic to form bubbles.

SHERAPRESS does not adhere to highly cross-linked artificial teeth. Artificial teeth must be pretreated (roughened) and/or provided with mechanical retention.

Insulation

Plaster surfaces are wetted with two thin layers of alginate insulation (e.g. SHERATEC), ensuring that the first layer is not yet dry when the second is applied.

The alginate insulation must be dry before further processing with acrylic to prevent any white discolourations on the base surfaces.

Mixing ratio 10:7

Recommended mixing ratio: 10 g powder with 7 g liquid. Free dosing is possible but a guarantee of the physical material properties is only given if the components are weighed out.

Mixing

Place the liquid in the mixing bowl (made of PE, PP or silicone) and quickly sprinkle over the powder. If free dosing, sprinkle over enough powder to completely bind the liquid in the mixing bowl. Then thoroughly mix together with a broad spatula (made of PE, PP or metal). Ensure that there are no bubbles.

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Do not use the material after the expiry date has passed. Do not store above 25°C. Keep out of direct sunlight. Always seal the container well after use. Keep away from heat/sparks/open flames/hot surfaces.Do not smoke.

8. Disposal

Dispose of in accordance with official regulations. Do not dispose of individual components in household waste or the sewage system. Completely polymerised plastic, however, can be disposed of in household waste.

9. Feedback

Please quote batch number and article number in all correspondence about the product. According to the EU Medical Device Directive, users/patients are obligated to report serious incidents related to a medical device to the manufacturer and the national competent authorities in the country in which they occurred.

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